

## history log for data analysis

### 1. converting DICOMS to BIDS data-structure

- get the scanner-generated and time-ordered DICOM-files relating to a measured subject (“time-folder”) from the server (USB-stick)
- create data folder (“raw\_data\_tnac”) and within folder subdirectory (sub...)
- use the heudiconf docker-image and follow the instructions on <http://nipy.org/heudiconv/> to transfer data into BIDS-structure (define a “unique heudiconf.py-file”; watch out, that subject folders are just named 01, 02,...; adjust paths to py-file, DICOM-images and output folder)

```
sudo docker run --rm -it -v /home/lmn/Desktop/rawdata_TNAC/:/data:ro -v /home/lmn/Desktop/TNAC_BIDS/:/output nipy/heudiconv:latest -d /data/{subject}/*/*/* -s 02 -f /data/TNAC_heuristic.py -b -o /output
```

(just adjust subject number “02” to use order on new subject-data, the heudiconv.py-file is only generated once and will be used for all subjects!)

```
sudo docker run --rm -it -v /media/lmn/86A406A0A406933B/rawdata_TNAC/:/data:ro -v /media/lmn/86A406A0A406933B/TNAC_BIDS/:/output nipy/heudiconv:latest -d /data/{subject}/*/*/* -s 02 -f /data/TNAC_heuristic.py -b -o /output
```

- use the online BIDS-validator to check, if your data is correctly formatted into BIDS-structure (<http://incf.github.io/bids-validator/>); use firefox and not opera as browser

### 2. defacing

- change read-write-execute rights of your BIDS-data (chmod... via bash)

```
sudo chmod -R 777 /home/lmn/Desktop/TNAC_BIDS/
```

```
sudo chmod -R 777 /media/lmn/86A406A0A406933B/TNAC_BIDS/
```

- use the pydeface.py module to deface your T1 and T2 images (both anatomical; T2\* is functional!; Set up a python2 environment in bash; “source activate...deactivate”)

```
source activate python2
```

```
pydeface.py Desktop/TNAC_BIDS/sub-01/anat/sub-01_T1w.nii.gz
```

(change “sub-x” and “T1 to T2” to deface both anatomicals in all subjects)

```
pydeface.py /media/lmn/86A406A0A406933B/TNAC_BIDS/sub-02/anat/sub-02_T1w.nii.gz
```

```
source deactivate
```

- control via freeview or fslview
- rename the defaced images so that they still fit into BIDS and remove the non-defaced images

### 3. quality control

- use the mriqc docker-image and follow the instructions on <http://mriqc.readthedocs.io/en/stable/docker.html> to perform quality control on your BIDS

```
sudo docker run -it --rm -v ~/Desktop/TNAC_BIDS:/data:ro -v  
~/Desktop/TNAC_BIDS/derivatives/mriqc:/out poldracklab/mriqc:latest /data /out  
--verbose-reports participant --no-sub --participant_label 01
```

(can be run in parallel for more participants: e.g. 01 02 03 at end of commandline)

```
sudo docker run -it --rm -v /media/lmn/86A406A0A406933B/TNAC_BIDS:/data:ro  
-v /media/lmn/86A406A0A406933B/TNAC_BIDS/derivatives/mriqc:/out  
poldracklab/mriqc:latest /data /out --verbose-reports participant --no-sub  
--participant_label 02 03
```

- visually inspect your images (html-output in created folder “reports” within derivatives → 2 anatomical T1, T2- and 4 functional run 1-4 files), more detailed explanations on <https://mriqc.readthedocs.io/en/latest/reports/smri.html> and <https://mriqc.readthedocs.io/en/latest/reports/bold.html>

### 4. anatomical “preprocessing” via FreeSurfer, ANTS, mindboggle

- use the mindboggle docker-image and follow the instructions on <http://mindboggle.readthedocs.io/en/latest/> to process your anatomical data

```
sudo docker run --rm -it -v /home/lmn/Desktop/TNAC_BIDS:/home/jovyan/work/data  
bids/mindboggle /home/jovyan/work/data /home/jovyan/work/data/derivatives/ participant  
--participant_label 02
```

- new image for 3 Tesla data!

```
HOST=/media/lmn/86A406A0A406933B/TNAC_BIDS/
```

```
DOCK=/home/jovyan/work
```

```
sudo docker run --rm -ti -v $HOST:DOCK --entrypoint /bin/bash  
bids/mindboggle
```

- define subjects in script and run:

```
bash TNAC_mindboggle.sh
```

- inspect output via fslview or freeview

### 5. “functional” preprocessing

- write preprocessing workflow in own jupyter notebook, not within docker images!
- (use as much SPM-interface-based nodes as possible for later DCM!!!)
- (segment to dissect grey matter from white matter, brainstem, cerebrospinal fluid (use templates from the freesurfer “anatomical processing”))

- unzip the functional T2-images and use the selectfiles-node for input-control
- (no slice time correction with short TR=1; no despike or TSNR, because we don't know influences on time series needed for DCM; no scan drop, because scanner already performed shim-correction)
- motion-correct the functional images via MCFLIRT to a mean func image and use the parameters-file to visually inspect for motion-outliers via ArtifactDetect
- skullstrip the mean func image (BET)
- coregister the skullstripped mean func image to the skullstripped anatomical image (from the freesurfer-derived anatomical processing); use bbregister from freesurfer to create the "coregistration-matrix"
- create a transformation-matrix/"normalisation-matrix" (registration between subject's structural and MNI template) with ANTS-Registration
- Convert the BBRegister-transformation (out\_fsl\_file) to ANTS ITK-format for later merging of coregistration matrix with transformation matrix from antsreg
- put the 2 matrices into a list and perform coregistration and normalisation in one step; use AppliedTransform via ANTs for normalisation (before 1<sup>st</sup> level, because DCM!); Apply Transformation - applies the normalization matrix to contrast images (functional)
- for more information about normalisation via ANTs see <http://miykael.github.io/nipype-beginner-s-guide/normalize.html> (complete transformation)
- smoothing with FWHM=6mm of the realigned and normalized functionals